Rat Nuclear receptor subfamily 4 group A member 2 (NR4A2) ELISA Kit

SAB Signalway Antibody

Catalog No: #EK8738

Package Size: #EK8738-1 48T #EK8738-2 96T

Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

ט	es	cr	pi	IIO	n

Rat Nuclear receptor subfamily 4 group A member 2 (NR4A2) ELISA Kit		
ELISA Kit		
ELISA		
Rat (Rattus norvegicus)		
HZF-3; NOT; NURR1; RNR1; TINUR; NGFI-B/nur77 beta-type transcription factor		
homolog OTTHUMP00000204707 T-cell nuclear receptor NOT intermediate-early receptor protein nur related		
protein-1; human ho		
Q07917		
Q07917		
54278;		
The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5%		
within the expiration date under appropriate storage condition.		
The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days,		
and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China		
Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage		
at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).		

Application Details

Detect Range:0.156-10 ng/mL		
Sensitivity:0.078 ng/mL		
Sample Type:Serum, Plasma, Other biological fluids		
Sample Volume: 1-200 μL		
Assay Time:1-4.5h		
Detection wavelength:450 nm		

Product Description

Detection Method:SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate NR4A2 in samples. An antibody specific for NR4A2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyNR4A2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for NR4A2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of NR4A2 bound in the initial step. The color development is stopped and the intensity of the color is measured.

Note: This product is for in vitro research use only