Human Anti-papillomavirus type 17 IgG antibody (HPV17-Ab-IgG) ELISA Kit

Catalog No: #EK11781

Description

Package Size: #EK11781-1 48T #EK11781-2 96T



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description				
Product Name	Human Anti-papillomavirus type 17 IgG antibody (HPV17-Ab-IgG) ELISA Kit			
Brief Description	ELISA Kit			
Applications	ELISA			
Species Reactivity	Human (Homo sapiens)			
Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5%			
	within the expiration date under appropriate storage condition.			
	The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days,			
	and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China			
	Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage			
	at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).			

Application Details			
Detect Range:Detection antibo	ly		
Sensitivity:Request Information			
Sample Type:Serum, Plasma,)ther biological fluids		
Sample Volume: 1-200 µL			
Assay Time:1-4.5h			
Detection wavelength:450 nm			

Product Description

Detection Method:Competitive ELISATest principle:This assay employs the competitive enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with an antibody specific to HPV17-Ab-IgG. Standards or samples are then added to the appropriate microtiter plate wells with a Horseradish Peroxidase (HRP)-conjugated HPV17-Ab-IgG and incubated. The competitive inhibition reaction is launched between with HRP labeled HPV17-Ab-IgG and unlabeled HPV17-Ab-IgG with the antibody. A substrate solution is added to the wells and the color develops in opposite to the amount of HPV17-Ab-IgG in the sample. The color development is stopped and the intensity of the color is measured.

Note: This product is for in vitro research use only