Mre11 (Phospho-Ser676) Conjugated Antibody

Catalog No: #C12731

Package Size: #C12731-AF350 100ul #C12731-AF405 100ul #C12731-AF488 100ul #C12731-AF555 100ul #C12731-AF594 100ul #C12731-AF647 100ul #C12731-AF680 100ul #C12731-AF750 100ul #C12731-Biotin 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description

| Product Name | Mre11 (Phospho-Ser676) Conjugated Antibody |
|-----------------------|---|
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Species Reactivity | Hu |
| Specificity | Mre11 (Phospho-Ser676) Antibody detects endogenous levels of Mre11 only when phosphorylated at Ser676 |
| Immunogen Description | A synthesized peptide derived from human Mre11 (Phospho-Ser676) |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | MRE11A, AT-like disease, Endo/exonuclease Mre11, HNGS1, MRE11, MRE11 homolog 1, MRE11 homolog |
| | A, MRE11B, ATLD |
| Accession No. | Swiss-Prot#:P49959NCBI Gene ID:4361NCBI mRNA#:NCBI Protein#: |
| Uniprot | P49959 |
| GenelD | 4361; |
| Excitation Emission | AF350: 346nm/442nm |
| | AF405: 401nm/421nm |
| | AF488: 493nm/519nm |
| | AF555: 555nm/565nm |
| | AF594: 591nm/614nm |
| | AF647: 651nm/667nm |
| | AF680: 679nm/702nm |
| | AF750: 749nm/775nm |
| Calculated MW | 80 |
| Formulation | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide |
| Storage | Store at 4°Cin dark for 6 months |
| | |

Application Details

Suggested Dilution: AF350 conjugated: most applications: 1: 50 - 1: 250 AF405 conjugated: most applications: 1: 50 - 1: 250 AF488 conjugated: most applications: 1: 50 - 1: 250 AF555 conjugated: most applications: 1: 50 - 1: 250 AF594 conjugated: most applications: 1: 50 - 1: 250 AF647 conjugated: most applications: 1: 50 - 1: 250 AF680 conjugated: most applications: 1: 50 - 1: 250 AF750 conjugated: most applications: 1: 50 - 1: 250

Product Description

Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatogramphy using non-phosphopeptide.

Note: This product is for in vitro research use only