#### **Product Datasheet**

# FLT3 (Phospho-Tyr842) Conjugated Antibody

Catalog No: #C11731



Package Size: #C11731-AF350 100ul #C11731-AF405 100ul #C11731-AF488 100ul #C11731-AF555 100ul #C11731-AF594 100ul #C11731-AF647 100ul #C11731-AF680 100ul #C11731-AF750 100ul #C11731-Biotin 100ul

Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

# Description

Product Name	FLT3 (Phospho-Tyr842) Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of FLT3 only when phosphorylated at tyrosine 842.
Immunogen Description	Peptide sequence around phosphorylation site of tyrosine 842(S-N-Y(p)-V-V) derived from Human FLT3 .
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD135;STK1;FLK-2;FLT-3
Accession No.	Swiss-Prot#:P36888NCBI Gene ID:2322NCBI mRNA#:NM_004119.2. NCBI Protein#:NP_004110.2.
Uniprot	P36888
GenelD	2322;
Excitation Emission	AF350: 346nm/442nm
	AF405: 401nm/421nm
	AF488: 493nm/519nm
	AF555: 555nm/565nm
	AF594: 591nm/614nm
	AF647: 651nm/667nm
	AF680: 679nm/702nm
	AF750: 749nm/775nm
Calculated MW	170
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

# **Application Details**

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250		
AF405 conjugated: most applications: 1: 50 - 1: 250		
AF488 conjugated: most applications: 1: 50 - 1: 250		
AF555 conjugated: most applications: 1: 50 - 1: 250		
AF594 conjugated: most applications: 1: 50 - 1: 250		
AF647 conjugated: most applications: 1: 50 - 1: 250		
AF680 conjugated: most applications: 1: 50 - 1: 250		
AF750 conjugated: most applications: 1: 50 - 1: 250		
Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000		

### **Product Description**

Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

### Background

FLT3 encodes a class III receptor tyrosine kinase that regulates hematopoiesis. The receptor consists of an extracellular domain composed of five immunoglobulin-like domains, one transmembrane region, and a cytoplasmic kinase domain split into two parts by a kinase-insert domain. The receptor is activated by binding of the fms-related tyrosine kinase 3 ligand to the extracellular domain, which induces homodimer formation in the plasma membrane leading to autophosphorylation of the receptor. The activated receptor kinase subsequently phosphorylates and activates multiple cytoplasmic effector molecules in pathways involved in apoptosis, proliferation, and differentiation of hematopoietic cells in bone marrow. Mutations that result in the constitutive activation of this receptor result in acute myeloid leukemia and acute lymphoblastic leukemia.

Note: This product is for in vitro research use only