NOS2 Conjugated Antibody

Catalog No: #C35843

SAB Signalway Antibody

Package Size: #C35843-AF350 100ul #C35843-AF405 100ul #C35843-AF488 100ul

#C35843-AF555 100ul #C35843-AF594 100ul #C35843-AF647 100ul

#C35843-AF680 100ul #C35843-AF750 100ul #C35843-Biotin 100ul

Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description

· · · · · · · · · · · · · · · · · · ·	
Product Name	NOS2 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NOS2 protein.
Immunogen Description	Fusion protein corresponding to residues near the N terminal of human nitric oxide synthase 2, inducible
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	NOS; INOS; NOS2A; HEP-NOS
Accession No.	Swiss-Prot#:P35228NCBI Gene ID:4843NCBI Protein#:BC130283
Uniprot	P35228
GeneID	4843;
Excitation Emission	AF350: 346nm/442nm
	AF405: 401nm/421nm
	AF488: 493nm/519nm
	AF555: 555nm/565nm
	AF594: 591nm/614nm
	AF647: 651nm/667nm
	AF680: 679nm/702nm
	AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
AF405 conjugated: most applications: 1: 50 - 1: 250
AF488 conjugated: most applications: 1: 50 - 1: 250
AF555 conjugated: most applications: 1: 50 - 1: 250
AF594 conjugated: most applications: 1: 50 - 1: 250
AF647 conjugated: most applications: 1: 50 - 1: 250
AF680 conjugated: most applications: 1: 50 - 1: 250
AF750 conjugated: most applications: 1: 50 - 1: 250

 $Biotin \ conjugated: working \ with \ enzyme-conjugated \ streptavidin, \ most \ applications: \ 1:50 - 1:1,000$

Background

Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. This gene encodes a nitric oxide synthase which is expressed in liver and is inducible by a combination of lipopolysaccharide and certain cytokines. Three related pseudogenes are located within the Smith-Magenis syndrome region on chromosome 17.

Note: This product is for in vitro research use only