Product Datasheet

Ferritin light chain Conjugated Monoclonal Antibody

Catalog No: #C42045



Package Size: #C42045-AF350 100ul #C42045-AF405 100ul #C42045-AF488 100ul #C42045-AF555 100ul #C42045-AF594 100ul #C42045-AF647 100ul #C42045-AF680 100ul #C42045-AF750 100ul #C42045-Biotin 100ul

Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description

Product Name	Ferritin light chain Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Species Reactivity	Ни
Specificity	specific for human Feritin light chain denatured and native forms
Immunogen Description	Recombinant Human Ferritin light chain
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	FTL,Ferritin L subunit
Accession No.	Swiss-Prot#:P02792
Uniprot	P02792
GenelD	2512;
Excitation Emission	AF350: 346nm/442nm
	AF405: 401nm/421nm
	AF488: 493nm/519nm
	AF555: 555nm/565nm
	AF594: 591nm/614nm
	AF647: 651nm/667nm
	AF680: 679nm/702nm
	AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:	
AF350 conjugated: most applications: 1: 50 - 1: 250	
AF405 conjugated: most applications: 1: 50 - 1: 250	
AF488 conjugated: most applications: 1: 50 - 1: 250	
AF555 conjugated: most applications: 1: 50 - 1: 250	
AF594 conjugated: most applications: 1: 50 - 1: 250	
AF647 conjugated: most applications: 1: 50 - 1: 250	
AF680 conjugated: most applications: 1: 50 - 1: 250	
AF750 conjugated: most applications: 1: 50 - 1: 250	
Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000	

Background

Stores iron in a soluble, non-toxic, readily available form. Important for iron homeostasis. Iron is taken up in the ferrous form and deposited as ferric hydroxides after oxidation. Also plays a role in delivery of iron to cells. Mediates iron uptake in capsule cells of the developing kidney.

Note: This product is for in vitro research use only