Citrate synthetase Conjugated Antibody

Catalog No: #C49726

SAB Signalway Antibody

Package Size: #C49726-AF350 100ul #C49726-AF405 100ul #C49726-AF488 100ul

#C49726-AF555 100ul #C49726-AF594 100ul #C49726-AF647 100ul

#C49726-AF680 100ul #C49726-AF750 100ul #C49726-Biotin 100ul

Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description

•	
Product Name	Citrate synthetase Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Recombinant protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CISY_HUMAN antibody Citrate synthase antibody Citrate synthase, mitochondrial antibody citrate
	synthetase antibody Cs antibody EC 2.3.3 antibody EC 2.3.3.1 antibody
Accession No.	Swiss-Prot#:O75390
Uniprot	O75390
GeneID	1431;
Excitation Emission	AF350: 346nm/442nm
	AF405: 401nm/421nm
	AF488: 493nm/519nm
	AF555: 555nm/565nm
	AF594: 591nm/614nm
	AF647: 651nm/667nm
	AF680: 679nm/702nm
	AF750: 749nm/775nm
Calculated MW	52 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
AF405 conjugated: most applications: 1: 50 - 1: 250
AF488 conjugated: most applications: 1: 50 - 1: 250
AF555 conjugated: most applications: 1: 50 - 1: 250
AF594 conjugated: most applications: 1: 50 - 1: 250
AF647 conjugated: most applications: 1: 50 - 1: 250
AF680 conjugated: most applications: 1: 50 - 1: 250
AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Citrate synthase is found in nearly all cells capable of oxidative metabolism. The enzyme citrate synthase exists in nearly all living cells and stands as a pace-making enzyme in the first step of the citric acid cycle (or Krebs cycle). Citrate synthase is localized within eukaryotic cells in the mitochondrial matrix, but is encoded by nuclear DNA rather than mitochondrial. It is synthesized using cytoplasmic ribosomes, then transported into the mitochondrial matrix. Citrate synthase is commonly used as a quantitative enzyme marker for the presence of intact mitochondria.

Note: This product is for in vitro research use only