

Sheep Lactoferrin(LTF) ELISA Kit

Catalog No: #EK21696

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Sheep Lactoferrin(LTF) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Sheep |
| Storage | The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C). |

Application Details

Detect Range:3.12 ng/mL - 200 ng/mL

Sensitivity:1 ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 &micro;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

The microtiter plate provided in this kit has been pre-coated with Sheep LTF. Samples are then added to the appropriate microtiter plate wells with a Horseradish Peroxidase (HRP)-conjugated antibody preparation specific for Sheep LTF and incubated. Then substrate solutions are added to each well. Only those wells that contain Sheep LTF and enzyme-conjugated Sheep LTF will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm.

Note: This product is for in vitro research use only