## Factor X Antibody FITC Conjugated

Catalog No: #C08006F



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description	Support: tech@signalwayantibody.co
Product Name	Factor X Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Activated factor Xa heavy chain
Conjugates	FITC
Target Name	Factor X
Other Names	Activated factor Xa heavy chain; Coagulation factor; Coagulation factor X; EC 3.4.21.6; F10 antibody
	FA10_HUMAN; Factor Xa; FX; FXA; OTTHUMP00000018735; Prothrombinase; Stuart factor; Stuart Prower
	factor; Stuart-Prower factor; Coagulation factor X; Factor X heavy chain.
Accession No.	NCBI Gene ID2159
Uniprot	P00742
GeneID	2159;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## **Application Details**

IF=1:50-200

## Background

Hemostasis following tissue injury involves the deployment of essential plasma procoagulants (Prothrombin and Factors X, IX, V and VIII), which are involved in a blood coagulation cascade that leads to the formation of insoluble Fibrin clots and the promotion of platelet aggregation. Coagulation Factor X (Stuart Prower factor, FX, F10) is a vitamin K-dependent, single chain serine protease that is synthesized in the liver and circulates as an inactive precursor. The mature form of Factor X (Factor X A) is generated by Factor IX A- or Factor VII A-mediated cleavage at the tripeptide sequence, Arg-Lys-Arg, to yield a disulfide linked dimer. Together with the cofactor Factor V A and Ca2+ on the surface of platelets or endothelial cells, Factor X A coordinates as part of the prothrombinase complex, which mediates proteolysis of Prothrombin into active Thrombin. Mutations at the Factor X locus resulting in Factor X deficiencies can contribute to hemorrhagic diathesis.

Note: This product is for in vitro research use only