

ADAMTS6 Antibody

Catalog No: #37083

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ADAMTS6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ADAMTS6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human ADAM metalloproteinase with thrombospondin type 1 motif, 6
Target Name	ADAMTS6
Other Names	ADAM-TS6; ADAMTS-6; ADAM-TS 6
Accession No.	Swiss-Prot#: Q9UKP5NCBI Gene ID: 11174Gene Accssion: NP_922932
Uniprot	Q9UKP5
GeneID	11174;
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:15-1:50

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #37083 at dilution 1/20.

Background

This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motifs) protein family. Members of the family share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin

type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains.

Note: This product is for in vitro research use only