

ADHFE1 Antibody

Catalog No: #37085

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ADHFE1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total ADHFE1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human alcohol dehydrogenase, iron containing, 1
Target Name	ADHFE1
Other Names	HOT; ADH8; HMFT2263
Accession No.	Swiss-Prot#: Q8IWW8NCBI Gene ID: 137872Gene Accssion: NP_653251
Uniprot	Q8IWW8
GeneID	137872;
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #37085 at dilution 1/50.

Background

The ADHFE1 gene encodes hydroxyacid-oxoacid transhydrogenase (EC 1.1.99.24), which is responsible for the oxidation of 4-hydroxybutyrate in mammalian tissues.

Note: This product is for in vitro research use only