

PDCD7 Antibody

Catalog No: #37209

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PDCD7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total PDCD7 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human Programmed cell death 7
Target Name	PDCD7
Other Names	ES18; HES18
Accession No.	Swiss-Prot#: Q8N8D1NCBI Gene ID: 10081Gene Accssion: NP_005698
Uniprot	Q8N8D1
GeneID	10081;
Concentration	2.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human colon cancer tissue using #37209 at dilution 1/60.



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #37209 at dilution 1/60.

Background

This gene encodes a 59 kDa protein that is associated with the U11 small nuclear ribonucleoprotein (snRNP), which is a component of the minor U12-type spliceosome responsible for catalyzing pre-mRNA splicing of U12-type introns.

Note: This product is for in vitro research use only