

TEP1 Antibody

Catalog No: #37271

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TEP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TEP1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the N terminal of human telomerase-associated protein 1
Target Name	TEP1
Other Names	TP1; TLP1; p240; TROVE1; VAULT2
Accession No.	Swiss-Prot#: Q99973NCBI Gene ID: 7011Gene Accssion: NP_009041
Uniprot	Q99973
GeneID	7011;
Concentration	1.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human ovarian cancer tissue using #37271 at dilution 1/50.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #37271 at dilution 1/50.

Background

This gene product is a component of the ribonucleoprotein complex responsible for telomerase activity which catalyzes the addition of new telomeres on the chromosome ends. The telomerase-associated proteins are conserved from ciliates to humans.

Note: This product is for in vitro research use only