

CNTNAP3 Antibody

Catalog No: #37498

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CNTNAP3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CNTNAP3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human contactin associated protein-like 3
Target Name	CNTNAP3
Other Names	CASPR3; CNTNAP3A; RP11-290L7.1; RP11-138L21.1
Accession No.	Swiss-Prot#: Q9BZ76NCBI Gene ID: 79937Gene Accssion: NP_387504
Uniprot	Q9BZ76
GeneID	79937;
Concentration	1.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human lymphoma tissue using #37498 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #37498 at dilution 1/40.

Background

The protein encoded by this gene belongs to the NCP family of cell-recognition molecules. This family represents a distinct subgroup of the neurexins. NCP proteins mediate neuron-glia interactions in vertebrates and glial-glia contact in invertebrates. The protein encoded by this gene may play a role in cell recognition within the nervous system. Alternatively spliced transcript variants encoding different isoforms have been described but their biological nature has not been determined.

Note: This product is for in vitro research use only