

DAB2IP Antibody

Catalog No: #37520

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DAB2IP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total DAB2IP protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human DAB2 interacting protein
Target Name	DAB2IP
Other Names	AIP1; AIP-1; AF9Q34; DIP1/2
Accession No.	Swiss-Prot#: Q5VWQ8NCBI Gene ID: 153090Gene Accssion: XP_005251776
Uniprot	Q5VWQ8
GeneID	153090;
Concentration	2.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using #37520 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #37520 at dilution 1/40.

Background

DAB2IP is a Ras (MIM 190020) GTPase-activating protein (GAP) that acts as a tumor suppressor. The DAB2IP gene is inactivated by methylation in prostate and breast cancers

Note: This product is for in vitro research use only