

CLEC6A Antibody

Catalog No: #37525

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLEC6A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CLEC6A protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human C-type lectin domain family 6, member A
Target Name	CLEC6A
Other Names	CLEC4N; CLECSF10
Accession No.	Swiss-Prot#: Q6EIG7NCBI Gene ID: 93978Gene Accssion: NP_001007034
Uniprot	Q6EIG7
GeneID	93978;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #37525 at dilution 1/20.

Background

Like dectin-1 (CLEC7A; MIM 606264), CLEC6A, or dectin-2, is a type II membrane receptor with an extracellular C-type lectin-like domain fold. However, unlike dectin-1, dectin-2 lacks an immunoreceptor tyrosine-based activation motif (ITAM) in its cytoplasmic domain

Note: This product is for in vitro research use only