

EPPIN Antibody

Catalog No: #37552

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	EPPIN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total EPPIN protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human epididymal peptidase inhibitor
Target Name	EPPIN
Other Names	CT71; CT72; WAP7; WFDC7; SPINLW1; dJ461P17.2
Accession No.	Swiss-Prot#: O95925NCBI Gene ID: 57119Gene Accssion: NP_065131
Uniprot	O95925
GeneID	100526773;57119;
SDS-PAGE MW	15kd
Concentration	2.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Immunohistochemistry: 1:25-1:100

Images



Gel: 10%SDS-PAGE


Lysates (from left to right): NIH/3T3 cell, human testis and mouse fat tissue

Amount of lysate: 40ug per lane

Primary antibody: 1/400 dilution

Secondary antibody dilution: 1/8000

Exposure time: 20 seconds



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #37552 at dilution 1/25.

Background

This gene encodes an epididymal protease inhibitor, which contains both kunitz-type and WAP-type four-disulfide core (WFDC) protease inhibitor consensus sequences. Most WFDC genes are localized to chromosome 20q12-q13 in two clusters: centromeric and telomeric. This gene is a member of the WFDC gene family and belongs to the telomeric cluster. The protein can inhibit human sperm motility, and polymorphisms in this gene are associated with male infertility. Read-through transcription also exists between this gene and the downstream WFDC6 (WAP four-disulfide core domain 6) gene.

Note: This product is for in vitro research use only