

GREB1 Antibody

Catalog No: #37610

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GREB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GREB1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human growth regulation by estrogen in breast cancer 1
Target Name	GREB1
Other Names	KIAA0575; Protein GREB1
Accession No.	Swiss-Prot#: Q4ZG55NCBI Gene ID: 9687Gene Accssion: NP_055483
Uniprot	Q4ZG55
GeneID	9687;
Concentration	2.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human colon cancer tissue using #37610 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #37610 at dilution 1/40.

Background

This gene is an estrogen-responsive gene that is an early response gene in the estrogen receptor-regulated pathway. It is thought to play an important role in hormone-responsive tissues and cancer. Three alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Note: This product is for in vitro research use only