

IL17D Antibody

Catalog No: #37647

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL17D Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total IL17D protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human interleukin 17D
Target Name	IL17D
Other Names	IL27; IL-27; IL-17D
Accession No.	Swiss-Prot#: Q8TAD2NCBI Gene ID: 53342Gene Accssion: NP_612141
Uniprot	Q8TAD2
GeneID	53342;
Concentration	3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human tonsil cancer tissue using #37647 at dilution 1/50.



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using #37647 at dilution 1/50.

Background

The protein encoded by this gene is a cytokine that shares the sequence similarity with IL17. The treatment of endothelial cells with this cytokine has been shown to stimulate the production of other cytokines including IL6, IL8 and CSF2/ GM-CSF. The increased expression of IL8 induced by this cytokine was found to be NF-kappa B-dependent.

Note: This product is for in vitro research use only