

PLCL1 Antibody

Catalog No: #37831

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PLCL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PLCL1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the N terminal of human phospholipase C-like 1
Target Name	PLCL1
Other Names	PLCE; PLCL; PRIP; PLDL1; PPP1R127
Accession No.	Swiss-Prot#: Q15111NCBI Gene ID: 5334Gene Accssion: NP_006217/Q15111
Uniprot	Q15111
GeneID	5334;
Concentration	2.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human prostate cancer tissue using #37831 at dilution 1/40.

Background

PLC-L (for PLC-deleted in lung carcinoma) is a PLC family member, the expression of which is decreased in small cell lung carcinomas. PLC-L is normally expressed in a variety of fetal and adult organs including the lung. PLC-L is thought to be involved in inositol phospholipid-based intracellular signaling cascade and aberrant expression of PLC-L may contribute to the genesis or progression of human lung carcinoma.

Note: This product is for in vitro research use only