

RIPK4 Antibody

Catalog No: #37870

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RIPK4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RIPK4 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human receptor-interacting serine-threonine kinase 4
Target Name	RIPK4
Other Names	DIK; PKK; PPS2; RIP4; ANKK2; NKRD3; ANKRD3
Accession No.	Swiss-Prot#: P57078NCBI Gene ID: 54101Gene Accssion: NP_065690/P57078
Uniprot	P57078
GeneID	54101;
Concentration	2.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #37870 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #37870 at dilution 1/40.

Background

The protein encoded by this gene is a serine/threonine protein kinase that interacts with protein kinase C-delta. The encoded protein can also activate NFkappaB and is required for keratinocyte differentiation. This kinase undergoes autophosphorylation.

Note: This product is for in vitro research use only