

HSV Mouse Monoclonal Antibody

Catalog No: #38056

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HSV Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification using immunogen.
Applications	WB
Specificity	The HSV?tag antibody can recognize?HSV tag fusion?proteins.
Target Name	HSV
Other Names	HSRG1; HSV-1 stimulation-related gene 1 protein; HSV-I stimulating-related protein; KIAA0872; MON1 homolog B (yeast)
Accession No.	Swiss-Prot#:Q7L1V2
Uniprot	Q7L1V2
GeneID	22879;
SDS-PAGE MW	60kd
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:5000~1:10000

Images



Western blot analysis of Recombinant protein, using #38056 diluted at 1) 1:5,000 2) 1:10,000

Background

HSV-tag peptide QPELAPEDPED conjugated to KLH. HSV- Tag antibody can recognize C-terminal, internal, and N-terminal HSV-tagged proteins.

Note: This product is for in vitro research use only