

LDHD Rabbit Polyclonal Antibody

Catalog No: #38105

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LDHD Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Rt Ms
Specificity	The LDHD Rabbit Polyclonal Antibody detects endogenous LDHD protein.
Target Name	LDHD
Other Names	D Lactate Dehydrogenase antibody, DLD antibody, Lactate dehydrogenase D antibody, MGC57726 antibody
Accession No.	Swiss-Prot#:Q86WU2
Uniprot	Q86WU2
GeneID	197257;
SDS-PAGE MW	55kd
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000~1:5000

Images



Western blot analysis of 1) HeLa, 2) Mouse Skeletal Muscle , 3) Mouse Liver, 4) Rat Kidney, using #38105 diluted at 1:5,000.

Background

LDHD belongs to the D-isomer specific 2-hydroxyacid dehydrogenase family. D-lactate and D-glycerate dehydrogenase activities have been found in a similar protein in yeast.

Note: This product is for in vitro research use only