

## AmCyan Mouse Polyclonal Antibody

Catalog No: #38080

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	AmCyan Mouse Polyclonal Antibody
Host Species	Mouse
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Rt Ms Mk Dg Chk Hm Rb Pg Sh
Specificity	The AmCyan Mouse Polyclonal Antibody detects AmCyan & AmCyan Tag protein.
Target Name	AmCyan
Other Names	AmCyan
Concentration	1.0mg/ml
Formulation	Mouse IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:2000~1:5000

## Images



Western blot analysis of AmCyan recombinant protein, using #38080 diluted at 1:5,000.

## Background

AmCyan has been adapted from the corresponding full length cDNA for higher solubility, brighter emission, and rapid chromophore maturation (8&uuml;C12 hours). AmCyan1 has been human-codon optimized for enhanced translation in mammalian cells. AmCyan is a very good blue fluorescent reporter protein, and useful in two-color analyses with DsRed2, DsRed-Express, or AsRed2, and in three-color analyses with ZsYellow and HcRed.

---

Note: This product is for in vitro research use only