

## EBFP Rabbit Polyclonal Antibody

Catalog No: #38087

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	EBFP Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Specificity	The EBFP Rabbit Polyclonal Antibody detects EBFP & EBFP Tag protein.
Immunogen Description	Recombinant protein corresponding to EBFP-tag.
Target Name	EBFP
Other Names	BFP; EBFP; Blue Fluorescent Protein; Enhanced Blue Fluorescent Protein
Concentration	1.0mg/ml
Formulation	IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:2000~1:5000

## Images



Western blot analysis of EBFP recombinant protein, using #38087 diluted at 1:5,000.

## Background

Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein “tags” to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

Note: This product is for in vitro research use only