

BHMT antibody

Catalog No: #22304

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	BHMT antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 213 of Human BHMT
Target Name	BHMT
Accession No.	Swiss-Prot:Q93088Gene ID:635
Uniprot	Q93088
GeneID	635;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 45kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug whole cell lysate)
 A: Hep G2
 10% SDS PAGE
 Primary antibody diluted at 1: 1000

Background

This gene encodes a cytosolic enzyme that catalyzes the conversion of betaine and homocysteine to dimethylglycine and methionine, respectively. Defects in this gene could lead to hyperhomocyst(e)inemia, but such a defect has not yet been observed. [provided by RefSeq]

Note: This product is for in vitro research use only