

BST2 antibody

Catalog No: #38318

Package Size: #38318-1 50ul #38318-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	BST2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total BST2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human BST2 (NP_004326.1).
Target Name	BST2
Other Names	BST2;CD317;TETHERIN
Accession No.	Uniprot:Q10589GeneID:684
Uniprot	Q10589
GeneID	684
SDS-PAGE MW	36kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using BST2 antibody.

Background

Bone marrow stromal cells are involved in the growth and development of B-cells. The specific function of the protein encoded by the bone marrow stromal cell antigen 2 is undetermined; however, this protein may play a role in pre-B-cell growth and in rheumatoid arthritis.

Note: This product is for in vitro research use only