

## GNA15 antibody

Catalog No: #38352

Package Size: #38352-1 50ul #38352-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GNA15 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GNA15 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human GNA15 (NP_002059.3).
Target Name	GNA15
Other Names	GNA15;GNA16
Accession No.	Uniprot:P30679GenelD:2769
Uniprot	P30679
GenelD	2769
SDS-PAGE MW	43KDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using GNA15 antibody.

## Background

---

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems.

---

Note: This product is for in vitro research use only