

FOXP1 antibody

Catalog No: #38637

Package Size: #38637-1 50ul #38637-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FOXP1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total FOXP1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human FOXP1 (NP_001012523.1).
Target Name	FOXP1
Other Names	FOXP1;12CC4;HSPC215;MFH;QRF1;hFKH1B
Accession No.	Uniprot:Q9H334GeneID:27086
Uniprot	Q9H334
GeneID	27086
SDS-PAGE MW	110kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using FOXP1 antibody.

Background

This gene belongs to subfamily P of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in the regulation of tissue- and cell type-specific gene transcription during both development and adulthood. Forkhead box P1 protein contains both DNA-binding- and protein-protein binding-domains. This gene may act as a tumor suppressor as it is lost in several tumor types and maps to a chromosomal region (3p14.1) reported to contain a tumor suppressor gene(s). Alternative splicing results in multiple transcript variants encoding different isoforms.

Note: This product is for in vitro research use only