

RAD54L2 antibody

Catalog No: #38733



Package Size: #38733-1 50ul #38733-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RAD54L2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total RAD54L2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human RAD54L2.
Target Name	RAD54L2
Other Names	ARIP4; HSPC325; SRISNF2L;
Accession No.	Swiss-Prot#: Q9Y4B4NCBI Gene ID: 23132
Uniprot	Q9Y4B4
GeneID	23132;
SDS-PAGE MW	162kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RAD54L2 at 1:1000 dilution.

Background

DNA helicase that modulates androgen receptor (AR)-dependent transactivation in a promoter-dependent manner. Not able to remodel mononucleosomes in vitro (By similarity).

Note: This product is for in vitro research use only