

GAB1 antibody

Catalog No: #38774

Package Size: #38774-1 50ul #38774-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GAB1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GAB1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human GAB1 (NP_997006.1).
Target Name	GAB1
Other Names	GAB1;GRB2
Accession No.	Uniprot:Q13480GeneID:2549
Uniprot	Q13480
GeneID	2549
SDS-PAGE MW	110kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of SH-SY5Y cells, using GAB1 antibody.

Background

The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only