

## BNIP2 antibody

Catalog No: #38791

Package Size: #38791-1 50ul #38791-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	BNIP2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total BNIP2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human BNIP2 (NP_004321.2).
Target Name	BNIP2
Other Names	BNIP2;BNIP-2;NIP2
Accession No.	Uniprot:Q12982GeneID:663
Uniprot	Q12982
GeneID	663
SDS-PAGE MW	45-55kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using BNIP2 antibody.

## Background

---

This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein, which protects cells from virally-induced cell death. The encoded protein also interacts with E1B 19 kDa-like sequences of BCL2, another apoptotic protector. Alternate splicing results in multiple transcript variants.

---

Note: This product is for in vitro research use only