

## FAIM3 antibody

Catalog No: #38823

Package Size: #38823-1 50ul #38823-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FAIM3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total FAIM3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human FAIM3.
Target Name	FAIM3
Other Names	FCMR; TOSO;
Accession No.	Swiss-Prot#: O60667NCBI Gene ID: 9214
Uniprot	O60667
GeneID	9214;
SDS-PAGE MW	43kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000IF □ 1:50 - 1:100

## Images



Western blot analysis of extracts of various cell lines, using FAIM3 at 1:1000 dilution.

## Background

---

Fc receptors specifically bind to the Fc region of immunoglobulins (Igs) to mediate the unique functions of each Ig class. FAIM3 encodes an Fc receptor for IgM (see MIM 147020) (Kubagawa et al., 2009)

---

Note: This product is for in vitro research use only