

CFDP1 antibody

Catalog No: #38826

Package Size: #38826-1 50ul #38826-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CFDP1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total CFDP1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human CFDP1.
Target Name	CFDP1
Other Names	p97; BCNT; CP27; SWC5; Yeti; CENP-29; BUCENTAUR;
Accession No.	Swiss-Prot#: Q9UEE9NCBI Gene ID: 10428
Uniprot	Q9UEE9
GeneID	10428;
SDS-PAGE MW	33kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CFDP1 at 1:1000 dilution.

Background

May play a role during embryogenesis.

Note: This product is for in vitro research use only