

## ZMYND11 antibody

Catalog No: #38827

Package Size: #38827-1 50ul #38827-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ZMYND11 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total ZMYND11 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human ZMYND11.
Target Name	ZMYND11
Other Names	BS69; BRAM1;
Accession No.	Swiss-Prot#: Q15326NCBI Gene ID: 10771
Uniprot	Q15326
GeneID	10771;
SDS-PAGE MW	70kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZMYND11 at 1:1000 dilution.

## Background

---

The protein encoded by this gene was first identified by its ability to bind the adenovirus E1A protein. The protein localizes to the nucleus. It functions as a transcriptional repressor, and expression of E1A inhibits this repression. Alternatively spliced transcript variants encoding different isoforms have been identified.

---

Note: This product is for in vitro research use only