

## INTS4 antibody

Catalog No: #38837

Package Size: #38837-1 50ul #38837-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	INTS4 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total INTS4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human INTS4 (NP_291025.3).
Target Name	INTS4
Other Names	INTS4;INT4;MST093
Accession No.	Uniprot:Q96HW7GenelD:92105
Uniprot	Q96HW7
GenelD	92105
SDS-PAGE MW	108KDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of Rat liver, using INTS4 antibody.



Western blot analysis of extracts of various cell lines, using INTS4 antibody.

## Background

INTS4 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005 [PubMed 16239144]).[supplied by OMIM, Mar 2008]

Note: This product is for in vitro research use only