

TCN1 antibody

Catalog No: #38893

Package Size: #38893-1 50ul #38893-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCN1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total TCN1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TCN1 (NP_001053.2).
Target Name	TCN1
Other Names	TCN1;HC;TC-1;TC1;TCI
Accession No.	Uniprot:P20061GenID:6947
Uniprot	P20061
GenID	6947
SDS-PAGE MW	48kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TCN1 antibody.

Background

This gene encodes a member of the vitamin B12-binding protein family. This family of proteins, alternatively referred to as R binders, is expressed in various tissues and secretions. This protein is a major constituent of secondary granules in neutrophils and facilitates the transport of cobalamin into cells.

Note: This product is for in vitro research use only