

AAAS antibody

Catalog No: #38905

Package Size: #38905-1 50ul #38905-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	AAAS antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total AAAS protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human AAAS (NP_056480.1).
Target Name	AAAS
Other Names	AAAS;AAA;AAASb;ADRACALA;ADRACALIN;ALADIN;GL003;aladin
Accession No.	Uniprot:Q9NRG9GenelD:8086
Uniprot	Q9NRG9
GenelD	8086
SDS-PAGE MW	59kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using AAAS antibody.

Background

The protein encoded by this gene is a member of the WD-repeat family of regulatory proteins and may be involved in normal development of the peripheral and central nervous system. The encoded protein is part of the nuclear pore complex and is anchored there by NDC1. Defects in this gene are a cause of achalasia-addisonianism-alacrima syndrome (AAAS), also called triple-A syndrome or Allgrove syndrome. Two transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only