

BLZF1 antibody

Catalog No: #38908

Package Size: #38908-1 50ul #38908-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BLZF1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total BLZF1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human BLZF1 (NP_003657.1).
Target Name	BLZF1
Other Names	BLZF1;GOLGIN-45;JEM-1;JEM-1s;JEM1;golgin-45
Accession No.	Uniprot:Q9H2G9GeneID:8548
Uniprot	Q9H2G9
GeneID	8548
SDS-PAGE MW	45kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details


WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using BLZF1 antibody.

Western blot analysis of extracts of various cell lines, using BLZF1 antibody.



Background

Required for normal Golgi structure and for protein transport from the endoplasmic reticulum (ER through the Golgi apparatus to the cell surface).

Note: This product is for in vitro research use only