

NCR1 antibody

Catalog No: #38916



Package Size: #38916-1 50ul #38916-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NCR1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total NCR1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human NCR1.
Target Name	NCR1
Other Names	LY94; CD335; NKP46; NK-p46;
Accession No.	Swiss-Prot#: O76036NCBI Gene ID: 9437
Uniprot	O76036
GeneID	9437;
SDS-PAGE MW	34kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Background

Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.

Note: This product is for in vitro research use only