

GDA antibody

Catalog No: #38919

Package Size: #38919-1 50ul #38919-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GDA antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GDA protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human GDA (NP_004284.1).
Target Name	GDA
Other Names	GDA;CYPIN;GUANASE;NEDASIN;GAH
Accession No.	Uniprot:Q9Y2T3GeneID:9615
Uniprot	Q9Y2T3
GeneID	9615
SDS-PAGE MW	51kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using GDA Antibody.

Background

This gene encodes an enzyme responsible for the hydrolytic deamination of guanine. Studies in rat ortholog suggest this gene plays a role in microtubule assembly. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only