

## SPAG5 antibody

Catalog No: #38929

Package Size: #38929-1 50ul #38929-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SPAG5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total SPAG5 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human SPAG5.
Target Name	SPAG5
Other Names	MAP126; DEEPEST; hMAP126;
Accession No.	Swiss-Prot#: Q96R06NCBI Gene ID: 10615
Uniprot	Q96R06
GeneID	10615;
SDS-PAGE MW	134kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

IF 1:20 - 1:50

## Images



Immunofluorescence analysis of HeLa cells using SPAG5 .  
Blue: DAPI for nuclear staining.

## Background

---

This gene encodes a protein associated with the mitotic spindle apparatus. The encoded protein may be involved in the functional and dynamic regulation of mitotic spindles.

---

Note: This product is for in vitro research use only