

## SAP30BP antibody

Catalog No: #38937

Package Size: #38937-1 50ul #38937-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	SAP30BP antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total SAP30BP protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human SAP30BP.
Target Name	SAP30BP
Other Names	HTRG; HTRP; HCNGP;
Accession No.	Swiss-Prot#: Q9UHR5NCBI Gene ID: 29115
Uniprot	Q9UHR5
GeneID	29115;
SDS-PAGE MW	33kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using SAP30BP at 1:1000 dilution.

## Background

---

SAP30BP (SAP30 binding protein) is a protein-coding gene. Diseases associated with SAP30BP include obesity, and schizophrenia.

---

Note: This product is for in vitro research use only