

TFPT antibody

Catalog No: #38938

Package Size: #38938-1 50ul #38938-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TFPT antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TFPT protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human TFPT.
Target Name	TFPT
Other Names	FB1; amida; INO80F;
Accession No.	Swiss-Prot#: P0C1Z6NCBI Gene ID: 29844
Uniprot	P0C1Z6
GeneID	29844;
SDS-PAGE MW	28kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000IF □ 1:50 - 1:100

Images



Immunofluorescence analysis of U2OS cells using TFPT .

Background

Appears to promote apoptosis in a p53/TP53-independent manner.

Putative regulatory component of the chromatin remodeling INO80 complex which is involved in transcriptional regulation, DNA replication and probably DNA repair.

Note: This product is for in vitro research use only