

## INTS5 antibody

Catalog No: #39059

Package Size: #39059-1 50ul #39059-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	INTS5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total INTS5 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human INTS5.
Target Name	INTS5
Other Names	INT5; KIAA1698;
Accession No.	Swiss-Prot#: Q6P9B9NCBI Gene ID: 80789
Uniprot	Q6P9B9
GeneID	80789;
SDS-PAGE MW	107kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using INTS5 at 1:1000 dilution.

## Background

---

INTS5 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005)

---

Note: This product is for in vitro research use only