

SNAP91 antibody

Catalog No: #39149

Package Size: #39149-1 50ul #39149-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SNAP91 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Mouse,Rat
Specificity	The antibody detects endogenous level of total SNAP91 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human SNAP91.
Target Name	SNAP91
Other Names	CALM; AP180;
Accession No.	Swiss-Prot#: O60641NCBI Gene ID: 9892
Uniprot	O60641
GeneID	9892;
SDS-PAGE MW	92kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SNAP91 at 1:1000 dilution.

Background

Adaptins are components of the adapter complexes which link clathrin to receptors in coated vesicles. Clathrin-associated protein complexes are believed to interact with the cytoplasmic tails of membrane proteins, leading to their selection and concentration. Binding of AP180 to clathrin triskelia induces their assembly into 60-70 nm coats (By similarity).

Note: This product is for in vitro research use only