

## LMO4 antibody

Catalog No: #39194

Package Size: #39194-1 50ul #39194-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	LMO4 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total LMO4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human LMO4 (NP_006760.1).
Target Name	LMO4
Other Names	LMO4
Accession No.	Uniprot:P61968GeneID:8543
Uniprot	P61968
GeneID	8543
SDS-PAGE MW	/
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

IF 1:50 - 1:100

## Images



Immunofluorescence analysis of MCF-7 cells using LMO4 antibody.

## Background

---

This gene encodes a cysteine-rich protein that contains two LIM domains but lacks a DNA-binding homeodomain. The encoded protein may play a role as a transcriptional regulator or as an oncogene.

---

Note: This product is for in vitro research use only