

POU6F1 antibody

Catalog No: #39195

Package Size: #39195-1 50ul #39195-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	POU6F1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total POU6F1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human POU6F1.
Target Name	POU6F1
Other Names	BRN5; MPOU; TCFB1;
Accession No.	Swiss-Prot#: Q14863NCBI Gene ID: 5463
Uniprot	Q14863
GeneID	5463;
SDS-PAGE MW	33kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

WB □ 1:500 - 1:2000IF □ 1:50 - 1:100

Images



Immunofluorescence analysis of A-549 cells using POU6F1 .



Western blot analysis of extracts of various cell lines, using POU6F1 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 90s.

Background

Transcription factor that binds preferentially to a variant of the octamer motif (5'-ATGATAAT-3').

Note: This product is for in vitro research use only