

TOM1-like protein 2 antibody

Catalog No: #22411

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	TOM1-like protein 2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 394 and 457 of TOM1-like protein 2
Target Name	TOM1-like protein 2
Accession No.	Swiss-Prot:Q6ZVM7Gene ID:146691
Uniprot	Q6ZVM7
GeneID	146691;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 50kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: JurKat

B: NT2D1

7.5% SDS PAGE

Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only