

Annexin A10 Antibody

Catalog No: #39230

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Annexin A10 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Protein A purified
Applications	ELISA WB
Species Reactivity	Hu Rt Ms
Specificity	Annexin A10 antibody detects endogenous levels of Annexin A10 protein. Not tested with other proteins
Immunogen Type	Peptide-KLH
Immunogen Description	Synthetic peptide conjugated to KLH
Target Name	Annexin A10
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA WB: 1:200

Images



Note: This product is for in vitro research use only